

LEAD-LINED UNIT DOSE CABINET




The Unit Dose Cabinet is designed for hot labs with limited space.

The upper left compartment holds two sharps containers to facilitate decay rotation. A small chute with removable shielded lid allows syringe disposal into the front container. A large rear port allows safe removal of the second decayed container. The upper right compart-

ment is for bulk storage of unit dose ammo boxes. The shielded door provides access without additional exposure.

Side-by-side middle drawers can be used for storage of syringes, sources, pigs, radioisotopes and other small items requiring lead shielding.

Phantoms and flood sources can be stored in the bottom drawer which runs the full width of the cabinet.



Cabinet drawing dimensions can be accessed at www.biodex.com.

SPECIFICATIONS:

Dimensions: 30.5" w x 24" depth x 36.5" h (77.5 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Drawer Dimensions:

Bottom Drawer: I.D.: 23.5" w x 21.8" depth x 4.5" h (59.7 x 55.4 x 11.4 cm)

Left & Right Drawers: I.D.: 9.38" w x 17.88" depth x 4.5" h (23.8 x 45.4 x 11.4 cm)

Upper Doors: I.D.: 11.5" w x 17.62" depth x 11" h (29.2 x 44.7 x 27.9 cm)

Doors and Drawers: Key-locked

Countertop: Stainless steel with 4" backsplash and .5" spillproof lip

Finish: Powder coat

Weight:

244-120: 1156 lb (524.3 kg)

244-121: 1614 lb (732 kg)

244-122: 2328 lb (1056 kg)

244-120 Cabinet, Unit Dose, .25" lead

244-121 Cabinet, Unit Dose, .5" lead

244-122 Cabinet, Unit Dose, 1" lead

Related:

039-338 Monoject Container, 4 qt., 10/pkg
Fits 039-325, 039-326, 039-350 and 039-330

039-387 Monoject Container, 8 qt., 10/pkg
Fits 039-325, 039-326 and 039-330